




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Owl D2 Wide-Gel- Electrophoresesystem Standard 40 Each 1.5mm

64,00 €
zzgl. MwSt &
Versand



Beschreibung

Owl D2 Wide-Gel-Electrophoresesystem Owl D2 Wide-Gel Electrophoresis System Standard, 40 teeth, 1.0mm thick tooth
Eliminate comb assembly with Thermo Scientific™ Owl™ D2 Wide-Gel Electrophoresis System Combs. These heavy-
duty, one-piece combs are for use with Owl™ D2 Wide-Gel Electrophoresis Systems. Bezeichnung: Standard Zur
Verwendung mit: Thermo Scientific™ Owl D2 Wide-Gel Electrophoresis Systems Artikelbeschreibung: Standard; 40
teeth; 1.5mm thick Anzahl der Zähne: 40 Dicke der Zähne (metrisch): 1.5mm Zur Verwendung mit (Geräte): Owltrade;
D2 Wide-Gel Electrophoresis Systems Typ: Combs

Zusätzliche Informationen

Art.-Nr.	LD-3703
Hersteller (Marke)	Thermo Scientific

