



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific P1-CST Multiple Gradient Gel Caster, 10 teeth Standard 10-tooth comb - 1.5mm

48,40 €
zzgl. MwSt &
Versand



Beschreibung

PI-CST Multiple Gradient Gel Caster, 10 teeth Cast one to five gels simultaneously with Fisherbrand™ PI-CST Multiple Gradient Gel Caster, one of the simplest devices to cast gradient gels. With an external gradient maker, gels can be cast with ease and great reproducibility The design and format of this unit will provide years of service in your lab
Bezeichnung: Standard 10-tooth comb Zur Verwendung mit: PI-CST System Artikelbeschreibung: 10 teeth; Tooth thickness, 1.5mm Anzahl der Zähne: 10 Dicke (metrisch): 1.5mm Zur Verwendung mit (Geräte): Owltrade; PI-CST System Typ: Combs

Zusätzliche Informationen

| | |
|--------------------|-------------------|
| Art.-Nr. | LD-4841 |
| Hersteller (Marke) | Thermo Scientific |

