



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific P2-CST Multiple Gradient Gel Caster Comb, 25 teeth Blank glass

17,40 €
zzgl. MwSt &
Versand



Beschreibung

P2-CST Multiple Gradient Gel Caster Comb, 25 teeth Cast one to five gels simultaneously with the Fisherbrand™ P2-CST Multiple Gradient Gel Caster, one of the simplest devices to cast gradient gels. With an external gradient maker, gels can be cast with ease and great reproducibility The design and format of this unit will provide years of service in your lab Bezeichnung: Blank glass Zur Verwendung mit: P10DS-1 Dual Gel Vertical System Artikelbeschreibung: Blank glass Dicke der Zähne (metrisch): – Dicke (angloamerikanisch): 0.18 in. Zur Verwendung mit (Geräte): Fisherbrandtrade; P2-CST System Typ: Blank Glass

Zusätzliche Informationen

Art.-Nr.	LD-4868
Hersteller (Marke)	Thermo Scientific

