



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific P2-CST Multiple Gradient Gel Caster Comb, 25 teeth, 1.5mm tooth thickness

48,40 €
zzgl. MwSt &
Versand



Beschreibung

P2-CST Multiple Gradient Gel Caster Comb, 25 teeth, 1.5mm tooth thickness Cast one to five gels simultaneously with the Fisherbrand™ P2-CST Multiple Gradient Gel Caster, one of the simplest devices to cast gradient gels. With an external gradient maker, gels can be cast with ease and great reproducibility The design and format of this unit will provide years of service in your lab Bezeichnung: Standard 25-tooth comb Zur Verwendung mit: P2-CST System Artikelbeschreibung: Comb, 25 teeth; Tooth thickness, 1.5mm Anzahl der Zähne: 25 Dicke der Zähne (metrisch): 1.5mm Zur Verwendung mit (Geräte): Fisherbrandtrade; P2-CST System Typ: Combs

Zusätzliche Informationen

Art.-Nr.	LD-4871
Hersteller (Marke)	Thermo Scientific

