



Thermo Scientific AC 200-S7 230V

Artikelbilder

Beschreibung

AC 200-S7 (230V/50Hz) HEATING BATH CIRCULATOR STEEL, AC200-S7, Amb+13 TO +200 C, 4-7 L BATH VOLUME, 2.0 kW HEATER, FORCE SUCTION PUMP 20 L/MIN REMOTE SENSOR PORT; USB; RAMP PROGRAM; MULTIFUNCTION PORT; 8 LANGUAGES; UL/CSA,

Zusätzliche Informationen

Art.-Nr.	LD-4276
Hersteller (Marke)	Thermo Scientific
VGKL Nummer	105336107

