



Thermo Scientific AC150-A10B 230V

Artikelbilder

Beschreibung

AC 150-A10B (230V/50Hz) OPEN REFRIGERATION BATH CIRCULATOR STEEL, AC150-A10B, -10 to +100 DEGREE C, 250 W COOLING CAPACITY AT 20C, 17-30 L BATH VOLUME, 2.0 kW HEATER, FORCE SUCTION PUMP 20 L/MIN REMOTE SENSOR PORT; 8 LANGUAGES;

Zusätzliche Informationen

Art.-Nr.	LD-4268
Hersteller (Marke)	Thermo Scientific
VGKL Nummer	105335410

