



Thermo Scientific SC100-A25B 230V

Artikelbilder

Beschreibung

SC 100-A25B (230V/50Hz) OPEN REFRIGERATION BATH CIRCULATOR STEEL, SC100-A25B, -25 to +100 DEGREE C, 500 W COOLING CAPACITY AT 20C, 13-21 L BATH VOLUME, 2.0 kW HEATER, PUMP 17 L/MIN USB;

Zusätzliche Informationen

Art.-Nr.	LD-4228
Hersteller (Marke)	Thermo Scientific
VGKL Nummer	105332425

