



Thermo Scientific SC150-A10B 230V

Artikelbilder

Beschreibung

SC 150-A10B (230V/50Hz) OPEN REFRIGERATION BATH CIRCULATOR STEEL, SC150-A10B, -10 to +100 DEGREE C, 250 W COOLING CAPACITY AT 20C, 17-30 L BATH VOLUME, 2.0 kW HEATER, PUMP 17 L/MIN USB;

Zusätzliche Informationen

Art.-Nr.	LD-4244
Hersteller (Marke)	Thermo Scientific
VGKL Nummer	105333410

