



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657




Amarell® Satz Dichte-Aräometer, Typ M100 060-190, DIN 12791/BS 718

660,72 € Ab

483,89 €

**zzgl. MwSt &
Versand**



Amarell
Precision-Thermometers and Hydrometers

- General Purpose Thermometers
- Precision Thermometers
- Precision Thermometers in Sets
- Industrial Thermometers
- Cylinderthermometers
- Dial Thermometers
- Quartz Glass Thermometers
- Thermometers for Calorimetry
- Thermometers for Meteorology
- Thermometers for Special Purposes
- Thermometers with Standard Ground Joint
- Thermometers According to DIN, ASTM, IP, BS, S.T.P.T.C. and AFNOR
- AMM-Therm[®]** - non-Wetting Low-Hazard Precision Thermometers
- AMM-Gallium[®]** - non-Wetting Low-Hazard Precision Thermometers

- General Purpose Hydrometers
- Hydrometers According to National and International Standards
- Hydrometers for Mineral Oil and Liquefied Gas Testing
- Hydrometers for Milk, Water, According to Baumé and for Sugar Solutions (Saccharimeters)
- Density Hydrometers and Alcoholometers
- Hydrometers for Special Applications
- Fittings for Hydrometric Measurements

PRECISION

Amarell GmbH & Co.KG Linderstraße 43 D-81892 Kreuzweithelm Germany www.amarell.de

Beschreibung

Satz Dichte-Aräometer, Typ M100 060-190, DIN 12791/BS 718, 0,600-2,000:0,002g/cm³, Genauigkeit + 1 Skalenteilungswert, 250mm lang, Bezugstemperatur 20°C, ohne Thermometer, Ablesung unten, 14 Instrumente im Etui, eichfähig

Zusätzliche Informationen

Art.-Nr.	LL-1068
Hersteller (Marke)	Amarell

