



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Germany
+49 (0)6221 / 8442-44
<https://www.neolab.de/en>

VAT No.: DE 143 450 657



Eppendorf DeepWell-Rotor incl. 2 DeepWell-Platten- Gehänge | E-1898

Product Images



Description

DeepWell-Rotor Typ A-2-DWP, Eppendorf

inkl. 2 DeepWell-Platten-Gehänge, max. Beladungshöhe 85 mm.
Max. 3.700 UpM

Additional Information

Autoclavable	No
With closure/ lid	No
For centrifuge model	Eppendorf, 5804 Eppendorf, 5804 R
Rotor is two-rowed	No
Rotor type	Ausschwingrotor für Deep Well-Platten
Rotor(s) is/are suitable for	Deep Well-Platten, 2x
Adapter for additional containers/ attachment of accessory parts available as accessory	No