





neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Germany
+49 (0)6221 / 8442-44
<https://www.neolab.de/en>

VAT No.: DE 143 450 657



neoLab Four-port connector for tubing inner Ø 2.25 mm, 10 pcs/pack | 2-0987

Product Images



Description

Crisscross connector made of polypropylene, autoclavable until 121°C.

The cone has a special form, that the pipes can't slide down.
suitable for high pressure applications.

10 pcs/bag

Additional Information

Material	polypropylene (PP)
For inner hose diameter 1	2,25 mm
Color	white
Shape	cross-shaped
Autoclavable	Yes
In a set	No
Miniature version	No
Hose connector type	cross-connector
Symmetry	für gleiche Schlauch-DM