




neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Germany
+49 (0)6221 / 8442-44
<https://www.neolab.de/en>

VAT No.: DE 143 450 657



neoLab Grid insert for thermo box small, 15 x 500 ml bottles | 1-7542

Product Images



Description

Grid inserts for laboratory glass bottles, exactly adjusted to the dimensions of the thermal box.

Grid height: 8 cm

Maximum bottle diameter:

500 ml = 90 mm

Accessory for thermal boxes (order number: 1-7335)

Additional Information

For small bottle/ vial type	SML
Alphanumeric coding	No
Material	polypropylene (PP)
Height	80 mm
Length	500 mm
Width / Breadth	300 mm
Color	black
Suited for	bottles, round
For containers with diameter up to	90 mm
Capacity	15
For containers with volume up to	500 ml
For small tube/ micro vessel type	SML
Wells (W x L)	3 x 5
Drill holes according to customer wishes	No
For bottle type	laboratory bottles, round