



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e



Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Heidolph PharMed® 2-stopper tube, ID: 0.2mm - SWS: 0.9mm

**€24.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

2-stopper tubing (40 cm) for Cassette small Application: Ideal for medicine, laboratory and research Properties: High fatigue bending strength Non-toxic, biocompatible Very low gas permeability Well suited for acids and alkalis Material: polypropylene-based thermoplastic elastomer with plasticizers, opaque beige Complies with standards: USP Class VI, GLP, Pharmacopoea and Europaea Temperature range: -51 to +135 °C Sterilization: autoclavable or sterilizable with ethylene oxide or radiation Restriction: detachment of additives possible Application with: Acids: good Alkalis: good Solvents: not suitable Pressure: good Vacuum: excellent Viscous media: good Sterile media: excellent

## Additional Information

---

No.	BH-1638
Manufacturer (Brand)	Heidolph
EAN	4058072778064
Gross weight	0.01
Transport temperature	Room temperature
Net weight	0.01000 kg
Country of origin	USA
Customs tariff number	39173200

