



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e



Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph Polyblock for 1 x 250 ml flask

€247.00
plus VAT &
Shipping

Product Images



Description

Chemical-resistant attachment segment for a 250 ml round bottom flask. The segment is made of aluminum and is suitable for heating temperatures up to 260 °C. Together with the StarFish base plate, the PolyBlocks are part of the basic elements of the entire StarFish Multi reaction system. The PolyBlocks are top segments that can be combined with each other for reaction vessels of different sizes (max. 5 different inserts per StarFish system).

Additional Information

No.	BH-1351
Manufacturer (Brand)	Heidolph
VGKL number	104150200
EAN	4058072775575
Gross weight	1.5
Transport temperature	Room temperature
Net weight	0.60000 kg
Country of origin	Great Britain
Packing width	0.20000 m
Packing height	0.11000 m
Packing depth	0.21000 m
Packaging volume	0.00500000 ccm
Customs tariff number	84199085

