




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph Polyblock for 3 x 25 mm test tubes

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Chemical-resistant attachment segment for max. 3 test tubes with 25 mm diameter. The segment is made of aluminum and is suitable for heating temperatures up to 260 °C. Together with the StarFish base plate, the PolyBlocks are part of the basic elements of the entire StarFish multi-reaction system. The PolyBlocks are top segments that can be combined with each other for reaction vessels of different sizes (max. 5 different inserts per StarFish system).

Additional Information

No.	BH-1352
Manufacturer (Brand)	Heidolph
VGKL number	104158210
EAN	4058072775582
Gross weight	1
Transport temperature	Room temperature
Net weight	0.30000 kg
Country of origin	Great Britain
Packing width	0.20000 m
Packing height	0.11000 m
Packing depth	0.21000 m
Packaging volume	0.00500000 ccm
Customs tariff number	84199085

