



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph Reax 2 incl. universal adapter -- EU cable

€1,668.00
plus VAT &
Shipping

Product Images



Description

Overhead shaker with quick-clamping technology for 2 vessels for different applications and a load weight of 1 kg The speed can be individually and continuously adjusted from 20 to 100 rpm A universal adapter for the use of sealed vessels with a height of 50 - 160 mm is included in the scope of delivery for you Select the optional adapter for mixing up to 20 test tubes simultaneously as an accessory

Additional Information

No.	BH-1740
Manufacturer (Brand)	Heidolph
VGKL number	524130020
EAN	4058072778903
Gross weight	9.22
Transport temperature	Room temperature
Net weight	6.76000 kg
Country of origin	Germany
Region of origin	Bavaria
Packing width	0.50000 m
Packing height	0.28000 m
Packing depth	0.66000 m
Packaging volume	0.09200000 cc
Customs tariff number	84798200

