



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Heidolph retort clamp

**€63.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Standard laboratory clamp for holding glassware on the standard ground joint for perfect vertical alignment. The four jaws - lined with cork on the inside and aluminum on the outside - can securely hold glassware with a diameter of 2 to 85 mm.

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | BH-1396          |
| Manufacturer (Brand)  | Heidolph         |
| VGKL number           | 104158175        |
| EAN                   | 4058072776022    |
| Gross weight          | 0.155            |
| Transport temperature | Room temperature |
| Net weight            | 0.15500 kg       |
| Country of origin     | Great Britain    |
| Customs tariff number | 84799070         |

