



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph Tygon® food hose, ID: 4.8mm - SWS: 1.6mm

€32.10
plus VAT &
Shipping

Product Images



Description

Suitable for single and multi-channel pump operation (not suitable for use with Cassette small) Application: Food industry, ideally suited for products with high fat content Features: Extremely resistant to chemicals, e.g. suitable for use with polar solvents Contains no plasticizers and oils Particularly long service life Transparent tubing for improved visual inspection Extremely flexible Material: thermoplastic Complies with standards: USP Class VI, FDA (21 CFR 177.2600) and GLP Temperature range: -78 to +71 °C Sterilization: autoclavable 30 min at 1 bar, sterilizable by radiation or ethylene oxide Application with: Acids: very good Alkalis: very good Solvents: good Pressure: good Vacuum: good Viscous media: good Sterile media: good

Additional Information

No.	BH-1710
Manufacturer (Brand)	Heidolph
EAN	4058072778620
Gross weight	0.05
Transport temperature	Room temperature
Net weight	0.05000 kg
Country of origin	USA
Customs tariff number	39173200

