



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph Tygon® standard extension hose, ID: 0.5mm - SWS: 0.9mm

€9.80
plus VAT &
Shipping

Product Images



Description

Extension tubing for Cassette small Application: General laboratory applications Properties: Non-toxic, non-oxidizing Good resistance to acids, alkalis and inorganic media Very low gas permeability, long service life Material: thermoplastic soft PVC, transparent Complies with standards: FDA (21 CFR 177.2601) and GLP Temperature range: -50 to +75 °C Sterilization: autoclavable at 120 °C, 1 bar at 30 min (takes on a milky color) or with ethylene oxide Limitation: plasticizers may dissolve Application with: Acids: good Alkalis: good Solvents: not suitable Pressure: good Vacuum: good Viscous media: excellent Sterile media: conditional

Additional Information

No.	BH-1690
Manufacturer (Brand)	Heidolph
EAN	4058072778446
Gross weight	0.004
Transport temperature	Room temperature
Net weight	0.00400 kg
Country of origin	USA
Customs tariff number	39173200

