



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph vacuum regulator Hei-VAC Control

€1,745.00
plus VAT &
Shipping

Product Images



Description

For precise vacuum control of the house vacuum or valve-controlled vacuum pumps. Entry-level model for digital vacuum control Easy installation on the support rod of the cooler or as a "stand-alone" device next to the evaporator Aeration button to prevent foaming Graphic display of the entire vacuum process Control valve already included RS 232 interface Multi-language menu navigation Compatible with all Hei-VAP rotary evaporators that do not have an integrated vacuum control Complete with integrated ceramic diaphragm vacuum sensor, control valve and aeration valve. With plug-in power supply, ready for operation, with manual.

Additional Information

No.	BH-2074
Manufacturer (Brand)	Heidolph
EAN	4058072784089
Gross weight	3
Transport temperature	Room temperature
Net weight	1.70000 kg
Country of origin	Germany
Region of origin	Baden-Württemberg
Packing width	0.50000 m
Packing height	0.28000 m
Packing depth	0.51000 m
Packaging volume	0.07100000 cc
Customs tariff number	90322000

