



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Julabo test tube insert for 90 microliter tubes 30 x 11/12 mm D., stainless steel

€152.00
plus VAT &
Shipping

Product Images



Description

The test tube insert 11/12 mm ø made of stainless steel is suitable for temperatures up to 100°C. A maximum of 90 microliter tubes can be placed. When tempering solvent-containing or acidic samples above 60°C, we recommend using stainless steel test tube inserts.

Additional Information

| | |
|-----------------------|------------------|
| No. | JN-0235 |
| Manufacturer (Brand) | Julabo |
| EAN | 4058072845292 |
| Transport temperature | Room temperature |
| Net weight | 0.80000 kg |
| Country of origin | Germany |
| Customs tariff number | 73269098 |

