



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN statistics printer for KERN scales with RS-233 data interface

€460.00
plus VAT &
Shipping

Product Images



Description

For printout of weighing values, recipes etc. or for GLP logging of weighing data with printout of date and time (not in connection with KERN counting scales or KERN moisture analyzers). Number generator with statistical functions like number and sum of all measured values, minimum and maximum value of the population, range (Max-Min). Mean value, standard deviations, available units: mg, g, kg, pcs. Language versions: D, GB, FR Thermal printing Dimensions housing W×D×H 165×140×50 mm Operation via mains adapter: 230 V / 50 Hz, standard EU/GB/USA,standard Speed: 0,75 lines/s Paper strip: 112 mm Option: paper rolls (10 pieces), operation via mains adapter: 230 V / 50 Hz, standard EU/GB/USA,standard Tip: Buy scale and printer together, then you get the interface cable for free! Not suitable for KERN PFB

Additional Information

No.	KP-5631
Manufacturer (Brand)	KERN
VGKL number	105835004
EAN	4058072151157
Transport temperature	Room temperature
Color	light gray dark gray
Manufacturer	CORE

