



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN Stable pit frame, stainless steel, for models with weighing plate size W×D×H 1500×1250×81

**€1,150.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Stable pit frame for the installation of the weighing bridge for barrier-free drive-up
- **Material:** stainless steel
- Dimensions W×D×H 1335×1585×80 mm
- **Note:** For verified scales, the weighing bridge must be fixed to the floor. Optionally by means of a drive-on ramp, a pair of foot plates or a pit frame
- not retrofittable

## Additional Information

---

No.	KP-0450
Manufacturer (Brand)	KERN
EAN	4058072125721
Transport temperature	Room temperature
autoclavable at 121°C	No
Material	Steel lacquered
sterile	No
Wide	0 mm
Height	0 mm
Depth	0 mm
Manufacturer	CORE

