



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**KERN platform steel, painted 230x230x110 mm:
Max 3 kg: e=1 g: d=100 mg:**

**Price on
request
plus VAT &
Shipping**



Description

- Stainless steel weighing plate, base steel painted
- 1 load cell, aluminum, silicone coated, IP65, OIML-R60 approval for verification, class III, 3000 e
- Level indicator and foot screws for precise leveling of the scale as standard, ensuring the most accurate weighing results

Additional Information

No.	KP-4602
Manufacturer (Brand)	KERN
EAN	4058072142087
Transport temperature	Room temperature
autoclavable at 121°C	No
Material	Steel lacquered
sterile	No
Wide	230 mm
Height	110 mm
Depth	230 mm
Manufacturer	CORE

