



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Weighing beam Max 1500 kg: d=0.5 kg

**Price on
request**

**plus VAT &
Shipping**

Product Images



Description

Flexible solution for large, voluminous or long weighing goods High mobility: thanks to battery operation (optional) and compact, flat design suitable for use at several locations Weigh beam: steel, painted, 4 silicone-coated aluminum load cells, dust and splash water protection IP67. Weighing beam as component also available without evaluation unit: KERN KFA^a; Stable handles for transport of the weighing beam UFA-L: One roller per weighing beam for comfortable transport of the scale Evaluation unit: Details see table base incl. UFA-S: Model with shorter weighing beams, ideal for weighing compact goods or animals in transport boxes Protective hood included in deliveryDid you know? Our floor scales are delivered in a robust wooden box. This protects the high quality weighing technology from environmental influences and stress during transport. KERN – always one idea ahead Shipping by forwarding agency. Please ask for dimensions, gross weight, shipping costs

Additional Information

No.	KP-5562
Manufacturer (Brand)	KERN
EAN	4058072150679
Transport temperature	Room temperature
autoclavable at 121°C	No
Material	Steel lacquered
sterile	No
Wide	170 mm
Height	85 mm
Depth	170 mm
Manufacturer	CORE

