



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN stainless steel measuring cells, 10000 kg CT 10000-3P2,

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

- Accuracy according to OIML R60 C3
- CE and RoHS compliant
- Dust and splash water protection IP67 (according to EN 60529), hermetically welded
- Nickel-plated steel
- Application: Mass as well as pressure force measurements under harsh environmental conditions
- Suitable for platform scales, hopper scales, scales embedded in the floor and other weighing equipment
- 4-wire connection
- Nominal sensitivity: 3 mV/V
- Note: EX version, 6-wire connection and accuracy class C4 or C5 on request

## Additional Information

---

No.	KP-5770
Manufacturer (Brand)	KERN
EAN	4058072152420
Transport temperature	Room temperature
Manufacturer	CORE

