



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN stainless steel measuring cells, 1000 kg CT 1000-3P2,

**Price on
request
plus VAT &
Shipping**

Product Images



Description

- Accuracy according to OIML R60 C3
- CE and RoHS compliant
- Dust and splash water protection IP67 (according to EN 60529), hermetically welded
- Nickel-plated steel
- Application: Mass as well as pressure force measurements under harsh environmental conditions
- Suitable for platform scales, hopper scales, scales embedded in the floor and other weighing equipment
- 4-wire connection
- Nominal sensitivity: 3 mV/V
- Note: EX version, 6-wire connection and accuracy class C4 or C5 on request

Additional Information

No.	KP-5771
Manufacturer (Brand)	KERN
EAN	4058072152437
Transport temperature	Room temperature
Manufacturer	CORE

