



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN load cell 2000 kg

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

- Accuracy class C3
- Dust and splash water protection IP67 (according to EN 60529, hermetically sealed)
- Field of application: for tension and compression force measurements
- Suitable for hanging scales, hopper scales and other weighing devices as well as force measurements and test benches
- 6-wire connection
- Nominal characteristic value: 2 mV/V

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | KP-3919          |
| Manufacturer (Brand)  | KERN             |
| EAN                   | 4058072135973    |
| Transport temperature | Room temperature |
| Manufacturer          | CORE             |

