



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN load cell 100 kg

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

- Accuracy according to OIML R60 C3
- Dust and splash water protection IP68/IP69K (according to EN 60529), hermetically welded
- Stainless steel
- Field of application: mass as well as pressure force measurements under harsh environmental conditions
- Suitable for platform scales, checkweighing
- Maximum platform size 10-50: 400×400 kg
- Maximum platform size 100-500 kg: 800×800 mm
- 4-wire connection (10-50 kg)
- 6-wire connection (100-500 kg)
- Nominal characteristic value: 2 mV/V
- Note: Version according to OIML R60 C4 or C5 on request

## Additional Information

---

No.	KP-2601
Manufacturer (Brand)	KERN
EAN	4058072131180
Transport temperature	Room temperature
Manufacturer	CORE

