



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN calibrator for regular adjustment of the sound level meter

€770.00
plus VAT &
Shipping

Product Images



Description

Calibrator for periodic adjustment of sound level meter, Class 1, and for testing linearity of sound level meters
Applicable standards: IEC60942:2003 Class 1, ANSI S1.40-1984, GB/T 15173-1994 Output frequency 1 kHz (+/- 0.5 %)
Output sound pressure, selectable 94 dB or 114 dB (± 0.3 dB) Distortion factor ≤ 2 % Stabilization time ≤ 10 s The calibrator is suitable for 1/2" as well as 1/4" microphones (adapter included) according to standard IEC 61094-4

Additional Information

No.	KP-0161
Manufacturer (Brand)	KERN
EAN	4058072124670
Transport temperature	Room temperature
Manufacturer	CORE

