



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN External measuring head, 5 MHz, Ø 12 mm, for hot test materials: Measuring range (steel) 1-225 mm

**€270.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- for hot test materials: Measuring range (steel)
- - 1-225 mm at normal temperatures
- - 4-100 mm at temperatures up to 300 °C

## Additional Information

---

No.	KP-0633
Manufacturer (Brand)	KERN
EAN	4058072123512
Transport temperature	Room temperature
Manufacturer	CORE

