



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN drag pointer 285,

€45.00
plus VAT &
Shipping

Product Images



Description

- The drag pointer consists of a plastic ring with a Teflon layer that glides well, making it easy to slip over the scale tube
- **Application:** Push the drag pointer to the zero point before measuring. When the scale is loaded, the drag pointer moves with it and stops at the maximum force per measuring operation
- Suitable for tensile tests up to 30 kg/300 N

Additional Information

| | |
|-----------------------|------------------|
| No. | KP-3404 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072134082 |
| Transport temperature | Room temperature |
| Manufacturer | CORE |

