



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN tare bowl with retaining bracket

€88.00
plus VAT &
Shipping

Product Images



Description

- made of stainless steel. Ideal for weighing loose parts, fruit, vegetables, etc.. Maximum load 30 kg. Removable tare bowl for easy and hygienic cleaning. Suitable for models KERN CH 15K20, HCB 20K10, HCN 20K50IP, HDB-N, HDB-XL

Additional Information

No.	KP-3437
Manufacturer (Brand)	KERN
EAN	4058072134273
Transport temperature	Room temperature
Manufacturer	CORE

