



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## KERN Manual test stand, one guide column Max 500 N &lt;-&gt; 230 mm

**€450.00**  
**plus VAT &**  
**Shipping**



## Description

---

Vertically and horizontally operable Precise measuring result High safety level for repeat measurements Large base plate with high mounting flexibility Suitable for all SAUTER force gauges up to 500 N (not included) Digital length measuring unit - Measuring range: 200 mm - Readability: 0.01 mm - Zero setting possible - Pre-length manually adjustable Spindle height from base plate: 297 mm Spindle stroke (stroke of one rotation): 3.1 mm

## Additional Information

---

No.	KP-5543
Manufacturer (Brand)	KERN
EAN	4058072150488
Transport temperature	Room temperature
Manufacturer	CORE

