



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## KERN Manual test stand, one guide column for TVP-L BOSCH Max 500 N

**€440.00**  
**plus VAT &**  
**Shipping**



## Description

---

Fast and consistent measurement High safety level for repeat measurements Wide range of applications and precise measurement results Slide design for distance measurements Large base plate with high mounting flexibility Suitable for all SAUTER force gauges up to 500 N (not included) Digital length measuring unit - Measuring range: 0.01 mm - Zero setting possible - Pre-length manually adjustable Maximum slide height above base plate: 318 mm Maximum slide stroke with one lever movement: 78 mm

## Additional Information

---

No.	KP-5556
Manufacturer (Brand)	KERN
EAN	4058072150617
Transport temperature	Room temperature
Manufacturer	CORE

