



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Set for density determination of liquids and solids with density $\rho \geq 1$

€380.00
plus VAT &
Shipping

Product Images



Description

- Density determination of liquids and solids with density ≤ 1
- Suitable for all KERN analytical balances
- **Scope of delivery:** beaker, thermometer, sample holder, adapter, balance weights, sinker (stainless steel, 20 g), universal immersion basket for floating and sinking solids. The wire of the immersion basket has a diameter of 0.5 mm and thus meets the requirements of DIN EN ISO 1183-1 (method for determining the density of non-foamed plastics)

Additional Information

No.	KP-5598
Manufacturer (Brand)	KERN
EAN	4058072150976
Transport temperature	Room temperature
Manufacturer	CORE

