



Macherey & Nagel AQUADUR 5-25/Ds/100

Product Images



Description

AQUADUR test strips Gradation: 5 > 10 > 15 > 20 > 25 °d Pack of 100 strips 6 x 95 mm

Additional Information

No.	ML-6305
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	244691201
EAN	4058072964931
Gross weight	0.0406
Manufacturer EAN	4046681004490
Storage temperature	Raumtemperatur
Net weight	0.04060 kg
S-phrases	S 22
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.04000 m
sterile	No
Packing height	0.12300 m
Packing depth	0.04000 m
Packaging volume	0.00019700 ccm
Customs tariff number	38229000

