



Macherey & Nagel NANO Aluminum

Product Images



Description

NANOCOLOR Aluminum Rectangular Cell Test Measuring range: 0.01-1.00 mg/L Al³⁺ Pg. à 250 determinations

Additional Information

No.	ML-6510
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	804691802
EAN	4058072966942
Gross weight	0.616
Storage temperature	Raumtemperatur
Net weight	0.61600 kg
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.16000 m
sterile	No
Packing height	0.13000 m
Packing depth	0.11000 m
Packaging volume	0.00228800 ccm
Customs tariff number	38229000

