



Macherey & Nagel NANO Ammonium

Product Images



Description

NANOCOLOR Ammonium Rectangular Cell Test Measuring range: 0.01-2.0 mg/L NH₄-N 0.01-2.5 mg/L NH₄⁺ Pg. à 100 determinations __UN 3316 Chemistry Test Kit 9 II 0.100 L ADR/GGVSEB M11 (E) ADR 3.3.1/251: LQ PAX 960 10 Kg CAO 960 10 Kg

Additional Information

No.	ML-6512
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	804691805
GHS symbol	GHS05 GHS07 GHS08
EAN	4058072966966
Gross weight	0.2387
H-phrases	H314. H361f
Storage temperature	Raumtemperatur
Net weight	0.23870 kg
P phrases	P201. P260sh. P280sh. P303. 361. 353. P305. 351. 338. P310
R-phrases	R 34
S-phrases	S 26-37/39-45
UN number	3316
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.16000 m
sterile	No
Packing height	0.13000 m
Packing depth	0.11000 m
Packaging volume	0.00228800 ccm
Customs tariff number	38229000
Signal word	Gefahr

