



Macherey & Nagel NANO Silica

Product Images



Description

NANOCOLOR Silica Rectangular Cell Test Measuring range: 0.01-5.00 mg/L Si 0.002-0.100 mg/L Si 0.02-10.0 mg/L SiO₂ 0.005-0.200 mg/L SiO₂ Pg. à 250 determinations __UN 3316 Chemistry Test Kit 9 II 0.200 L ADR/GGVSEB M11 (E) ADR 3.3.1/251: LQ PAX 960 10 Kg CAO 960 10 Kg

Additional Information

No.	ML-6528
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	804691848
GHS symbol	GHS05 GHS07
EAN	4058072967147
Gross weight	0.436
H-phrases	H318
Storage temperature	Raumtemperatur
Net weight	0.43600 kg
P phrases	P280sh. P305. 351. 338. P310
R-phrases	R 31-36/38-41
S-phrases	S 26-39-45
UN number	3316
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.16000 m
sterile	No
Packing height	0.13000 m
Packing depth	0.11000 m
Packaging volume	0.00228800 ccm
Customs tariff number	38229000
Signal word	Gefahr

