



Macherey & Nagel NANO Manganese

Product Images



Description

NANOCOLOR Manganese Rectangular Cell Test Measuring range: 0.01-10.0 mg/L Mn Pg. à 250 determinations __UN 3316 Chemistry Test Kit 9 II 0.200 L ADR/GGVSEB M11 (E) ADR 3.3.1/251: LQ PAX 960 10 Kg CAO 960 10 Kg

Additional Information

No.	ML-6534
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	804691860
GHS symbol	GHS02 GHS05 GHS07 GHS08 GHS09
EAN	4058072967208
Gross weight	0.421
H-phrases	H314. H317. H351
Storage temperature	Raumtemperatur
Net weight	0.42100 kg
P phrases	P201. P260sh. P280sh. P303. 361. 353. P305. 351. 338. P310
R-phrases	R 34-40-43
S-phrases	S 22-26-36/37/39-45-61
UN number	3316
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.16000 m
sterile	No
Packing height	0.13000 m
Packing depth	0.11000 m
Packaging volume	0.00228800 ccm
Customs tariff number	38229000
Signal word	Gefahr

