



Macherey & Nagel NANO Nitrate

Product Images



Description

NANOCOLOR Nitrate Rectangular Cell Test Measuring range: 0.9-30.0 mg/L NO₃-N 4-140 mg/L NO₃- Pg. à 100 determinations __UN
3316 Chemistry Test Kit 9 II 0.400 L ADR/GGVSEB M11 (E) ADR 3.3.1/251: LQ PAX 960 10 Kg CAO 960 10 Kg

Additional Information

No.	ML-6537
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
GHS symbol	GHS02 GHS05 GHS07
EAN	4058072967239
Gross weight	1.286
H-phrases	H314
Storage temperature	Raumtemperatur
Net weight	1.28600 kg
P phrases	P260sh. P280sh. P303. 361. 353. P305. 351. 338. P310
R-phrases	R 35
S-phrases	S 16-26-30-45
UN number	3316
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.14500 m
sterile	No
Packing height	0.13100 m
Packing depth	0.19500 m
Packaging volume	0.00370400 ccm
Customs tariff number	38229000
Signal word	Gefahr

