



Macherey & Nagel Nanocontrol Nitrite

Product Images



Description

NANOCONTROL Standard Nitrite 2x30 mL standard solution/ 30 mL addition solution 100+ sufficient for 15 determinations

Additional Information

No.	ML-6738
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	804692568
EAN	4058072969240
Gross weight	0.2033
Storage temperature	Raumtemperatur
Net weight	0.20330 kg
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.12000 m
sterile	No
Packing height	0.03800 m
Packing depth	0.10100 m
Packaging volume	0.00046100 ccm
Customs tariff number	38229000

