



Macherey & Nagel QUANTOFIX Calcium

Product Images



Description

QUANTOFIX Calcium test strips 6 x 95 mm Measuring range: 0-10-25-50-100 mg/L Ca²⁺ Pg. à 60 determinations incl. reagents and reaction vessel __UN 3316 Chemistry test kit 9 II 0.030 L ADR/GGVSEB M11 (E) ADR 3.3.1/251: LQ PAX 960 10 Kg CAO 960 10 Kg

Additional Information

No.	ML-6336
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	244691324
GHS symbol	GHS03 GHS05 GHS07
EAN	4058072965228
Gross weight	0.186
Manufacturer EAN	4046681004032
H-phrases	H314
Storage temperature	Raumtemperatur
Net weight	0.18600 kg
P phrases	P260sh. P280sh. P303. 361. 353. P305. 351. 338. P310
R-phrases	R 35
S-phrases	S 22-26-37/39-45
UN number	3316
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.14800 m
sterile	No
Packing height	0.04000 m
Packing depth	0.12200 m
Packaging volume	0.00072200 ccm
Customs tariff number	38229000
Signal word	Gefahr

