



Macherey & Nagel QUANTOFIX Potassium

Product Images



Description

QUANTOFIX Potassium test strips 6 x 95 mm Measuring range: 0-200-400-700-1000-1500 mg/L K+ Pg. à 100 determinations incl. reagents and reaction vessel

Additional Information

No.	ML-6328
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	244691316
EAN	4058072965143
Gross weight	0.1949
Manufacturer EAN	4046681004117
Storage temperature	Raumtemperatur
Net weight	0.19490 kg
S-phrases	S 22
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.14800 m
sterile	No
Packing height	0.04000 m
Packing depth	0.12200 m
Packaging volume	0.00072200 ccm
Customs tariff number	38229000

