



Macherey & Nagel QUANTOFIX Molybdenum

Product Images



Description

QUANTOFIX Molybdenum test strips 6 x 95 mm Measuring range: 0-5-20-50-100-250 mg/L Mo6+ Pg. à 100 determinations incl. reagents and reaction vessel __UN 3316 Chemistry test kit 9 II 0.040 kg ADR/GGVSEB M11 (E) ADR 3.3.1/251: LQ PAX 960 10 Kg CAO 960 10 Kg

Additional Information

No.	ML-6337
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	244691325
EAN	4058072965235
Gross weight	0.1816
Manufacturer EAN	4046681004025
Storage temperature	Raumtemperatur
Net weight	0.18160 kg
S-phrases	S 22-26
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.14800 m
sterile	No
Packing height	0.04000 m
Packing depth	0.12200 m
Packaging volume	0.00072200 ccm
Customs tariff number	38229000

