



Macherey & Nagel QUANTOFIX Nickel

Product Images



Description

QUANTOFIX Nickel test strips 6 x 95 mm Measuring range: 0-10-25-50-100- 250-500-1000 mg/L Ni²⁺ Pg. à 100 determinations

Additional Information

No.	ML-6318
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	244691305
EAN	4058072965068
Gross weight	0.0411
Manufacturer EAN	4046681004636
Storage temperature	Raumtemperatur
Net weight	0.04110 kg
S-phrases	S 22
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.04000 m
sterile	No
Packing height	0.12000 m
Packing depth	0.04000 m
Packaging volume	0.00019200 ccm
Customs tariff number	38229000

