



Macherey & Nagel QUANTOFIX Nitrite

Product Images



Description

QUANTOFIX Nitrite test strips 6 x 95 mm Measuring range: 0-1-5-10-20-40-80 mg/NO₂⁻ Pg. à 100 determinations

Additional Information

No.	ML-6323
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	244691311
EAN	4058072164263
Gross weight	0.0443
Manufacturer EAN	4046681002632
Storage temperature	Raumtemperatur
Net weight	0.04430 kg
S-phrases	S 22
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.04000 m
sterile	No
Packing height	0.12000 m
Packing depth	0.04000 m
Packaging volume	0.00019200 ccm
Customs tariff number	38229000

