



Macherey & Nagel QUANTOFIX Nitrite/pH

Product Images



Description

QUANTOFIX Nitrite/pH test strips 6 x 95 mm Measuring range: 0-1-5-10-20-40-80 mg/L NO₂- pH 6.0-6.4-6.7-7.0-7.3-7.6-7.9-8.2-8.4-8.6-8.8-9.0-9.3-9.6 Pg. à 100 determinations

Additional Information

No.	ML-6349
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	104691338
EAN	4058072178468
Gross weight	0.035
Manufacturer EAN	4046681003905
Storage temperature	Raumtemperatur
Net weight	0.03500 kg
S-phrases	S 22
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.04300 m
sterile	No
Packing height	0.04000 m
Packing depth	0.12300 m
Packaging volume	0.00021200 ccm
Customs tariff number	38229000

