



Macherey & Nagel QUANTOFIX Sulfate

Product Images



Description

QUANTOFIX Sulfate test strips 6 x 95 mm Measuring range: 400-> 800-> 1200- > 1600 mg/L SO₄²⁻ Pg. à 100 determinations

Additional Information

No.	ML-6341
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	244691329
EAN	4058072965266
Gross weight	0.0429
Manufacturer EAN	4046681003981
Storage temperature	Raumtemperatur
Net weight	0.04290 kg
S-phrases	S 22
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.04000 m
sterile	No
Packing height	0.12000 m
Packing depth	0.04000 m
Packaging volume	0.00019200 ccm
Customs tariff number	38229000

